

CACNG2 Polyclonal Antibody



Catalog Number: 161925

Product name

CACNG2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CACNG2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 204-323 of human CACNG2 (NP_006069.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:100

Validations



Immunofluorescence - CACNG2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CACNG2 antibody . Blue: DAPI for nuclear staining.