

# HLA-DPB1 Polyclonal Antibody



Catalog Number: 161935

## Product name

HLA-DPB1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to HLA-DPB1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-225 of human HLA-DPB1 (NP\_002112.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations




Western blot - HLA-DPB1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HLA-DPB1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunohistochemistry - HLA-DPB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lymphonodus using HLA-DPB1 antibody at dilution of 1:100 (40x lens).



# HLA-DPB1 Polyclonal Antibody



Catalog Number: 161935

---

Immunohistochemistry - HLA-DPB1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using HLA-DPB1 antibody at dilution of 1:100 (40x lens).