

# SERPINA12 Polyclonal Antibody



Catalog Number: 161948

---

## Product name

SERPINA12 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SERPINA12

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 20-320 of human SERPINA12 (NP\_776249.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



Immunofluorescence - SERPINA12 Polyclonal Antibody

Immunofluorescence analysis of PC-12 cells using SERPINA12 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.