

# UBIAD1 Polyclonal Antibody



Catalog Number: 161973

---

## Product name

UBIAD1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to UBIAD1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human UBIAD1 (NP\_037451.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - UBIAD1 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using UBIAD1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.