

ZNF157 Polyclonal Antibody



Catalog Number: 162045

Product name

ZNF157 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ZNF157

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 70-165 of human ZNF157 (NP_003437.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB

Dilutions

WB 1:500 - 1:2000

Validations

U-87MG



Western blot - ZNF157 Polyclonal Antibody

Western blot analysis of extracts of U-87MG cells, using ZNF157 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.