

CCDC59 Polyclonal Antibody



Catalog Number: 162068

Product name

CCDC59 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to CCDC59

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-241 of human CCDC59 (NP_054886.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - CCDC59 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using CCDC59 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.



Immunohistochemistry - CCDC59 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat liver using CCDC59 antibody at dilution of 1:100



(40x lens).


CCDC59 Polyclonal Antibody



Catalog Number: 162068

Immunohistochemistry - CCDC59 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat testis using CCDC59 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CCDC59 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using CCDC59 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - CCDC59 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human lung using CCDC59 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - CCDC59 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using CCDC59 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CCDC59 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon using CCDC59 antibody at dilution of




1:100 (40x lens).

Immunohistochemistry - CCDC59 Polyclonal Antibody

CCDC59 Polyclonal Antibody

Catalog Number: 162068


Immunohistochemistry of paraffin-embedded human uterine cancer using CCDC59 antibody at dilution of 1:100 (40x lens).



Immunofluorescence - CCDC59 Polyclonal Antibody

Immunofluorescence analysis of NIH-3T3 cells

using CCDC59 antibody at dilution of 1:100. Blue:



DAPI for nuclear staining.

Immunofluorescence - CCDC59 Polyclonal Antibody

Immunofluorescence analysis of U-2 OS cells using CCDC59 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.