

# ASCL4 Polyclonal Antibody



Catalog Number: 162087

## Product name

ASCL4 Polyclonal Antibody

## Specificity

Human, Rat

## Antibody description

Polyclonal antibody to ASCL4

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-172 of human ASCL4 (NP\_982260.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

## Validations



Western blot - ASCL4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ASCL4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.