

OTUD6B Polyclonal Antibody



Catalog Number: 162118

Product name

OTUD6B Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to OTUD6B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 31-323 of human OTUD6B (NP_057107.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IP

Dilutions

WB 1:500 - 1:2000


IP 1:50 - 1:200

Validations



Western blot - OTUD6B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using OTUD6B antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL



Immunoprecipitation - OTUD6B Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of U-87MG cells using 3ug OTUD6B antibody . Western blot was performed from the immunoprecipitate using OTUD6B antibody at a dilution of 1:1000.