

FLCN Polyclonal Antibody

Catalog Number: 162125

Product name

FLCN Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to FLCN

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-130 of human FLCN (NP_653207.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - FLCN Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using FLCN antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 10s.



Immunofluorescence - FLCN Polyclonal Antibody

Immunofluorescence analysis of C6 cells using FLCN antibody at dilution of 1:100. Blue: DAPI for




nuclear staining.

FLCN Polyclonal Antibody

Catalog Number: 162125

Immunofluorescence - FLCN Polyclonal Antibody

Immunofluorescence analysis of NIH/3T3 cells using FLCN antibody at dilution of 1:100. Blue:



DAPI for nuclear staining.

Immunofluorescence - FLCN Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using FLCN antibody at dilution of 1:100. Blue: DAPI for nuclear staining.