

ACHE Polyclonal Antibody



Catalog Number: 162135

Product name

ACHE Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to ACHE

Preparation

Antigen: A synthetic peptide of human ACHE

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF, FC

Dilutions


WB 1:200 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

FC 1:20 - 1:50

Validations



Western blot - ACHE Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using ACHE antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.