

# MYCT1 Polyclonal Antibody

Catalog Number: 162137

## Product name

MYCT1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to MYCT1

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 99-235 of human MYCT1 (NP\_079383.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations




## Western blot - MYCT1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using MYCT1 antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



## Immunofluorescence - MYCT1 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using MYCT1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

# MYCT1 Polyclonal Antibody



Catalog Number: 162137

---

Immunofluorescence - MYCT1 Polyclonal Antibody  
Immunofluorescence analysis of HeLa cells using

MYCT1 antibody at dilution of 1:100. Blue: DAPI  
for nuclear staining.