

HIGD1A Polyclonal Antibody



Catalog Number: 162154

Product name

HIGD1A Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to HIGD1A

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-93 of human HIGD1A (NP_054775.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

Validations



Western blot - HIGD1A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using HIGD1A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.