

MEAF6 Polyclonal Antibody



Catalog Number: 162184

Product name

Polyclonal

MEAF6 Polyclonal Antibody

Ig Type**Specificity**

Rabbit IgG

Human, Mouse

Applications**Antibody description**

WB

Polyclonal antibody to MEAF6

Dilutions**Preparation**

WB 1:500 - 1:2000

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-191 of human MEAF6 (NP_001257804.1).

Validations**Formulation**

Western blot - MEAF6 Polyclonal Antibody

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Western blot analysis of extracts of various cell lines, using MEAF6 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

Clonality