

IDH3A Polyclonal Antibody



Catalog Number: 162198

Product name

IDH3A Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to IDH3A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-366 of human IDH3A (NP_005521.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

Validations




Western blot - IDH3A Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using IDH3A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - IDH3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using IDH3A antibody at dilution of 1:100



(40x lens).


IDH3A Polyclonal Antibody



Catalog Number: 162198


Immunohistochemistry - IDH3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human thyroid cancer using IDH3A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - IDH3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using IDH3A antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - IDH3A Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using IDH3A antibody at dilution of 1:100 (40x lens).