

# ALOX12 Polyclonal Antibody



Catalog Number: 162224

## Product name

ALOX12 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ALOX12

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 483-663 of human ALOX12 (NP\_000688.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

## Validations



Immunohistochemistry - ALOX12 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using ALOX12 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ALOX12 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human colon carcinoma using ALOX12 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ALOX12 Polyclonal Antibody

# ALOX12 Polyclonal Antibody



Catalog Number: 162224

---

Immunohistochemistry of paraffin-embedded

mouse brain using ALOX12 antibody at dilution of 1:100 (40x lens).