

# DVL2 Polyclonal Antibody



Catalog Number: 162233

---

## Product name

DVL2 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DVL2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human DVL2 (NP\_004413.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


IHC, IF

## Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

## Validations



Immunofluorescence - DVL2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using DVL2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.