

ALDH6A1 Polyclonal Antibody



Catalog Number: 162243

Product name

ALDH6A1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ALDH6A1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 326-535 of human ALDH6A1 (NP_005580.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations



Immunohistochemistry - ALDH6A1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat ovary using ALDH6A1 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ALDH6A1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using ALDH6A1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ALDH6A1 Polyclonal Antibody

ALDH6A1 Polyclonal Antibody



Catalog Number: 162243

Immunohistochemistry of paraffin-embedded

mouse kidney using ALDH6A1 antibody at dilution of 1:100 (40x lens).