

NFE2L1 Polyclonal Antibody



Catalog Number: 162245

Product name

NFE2L1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to NFE2L1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 515-772 of human NFE2L1 (NP_003195.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Western blot - NFE2L1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using NFE2L1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.



Immunohistochemistry - NFE2L1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat spleen using NFE2L1 antibody at dilution of 1:100



(40x lens).


NFE2L1 Polyclonal Antibody



Catalog Number: 162245

Immunohistochemistry - NFE2L1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat heart using NFE2L1 antibody at dilution of 1:100



(40x lens).


Immunohistochemistry - NFE2L1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using NFE2L1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NFE2L1 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using NFE2L1 antibody at dilution of



1:100 (40x lens).


Immunohistochemistry - NFE2L1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human placenta using NFE2L1 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - NFE2L1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using NFE2L1 antibody at dilution of 1:100 (40x lens).




Immunohistochemistry - NFE2L1 Polyclonal Antibody

NFE2L1 Polyclonal Antibody

Catalog Number: 162245

Antibody

Immunohistochemistry of paraffin-embedded human uterus using NFE2L1 antibody at dilution




of 1:100 (40x lens).

Immunohistochemistry - NFE2L1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded

mouse brain using NFE2L1 antibody at dilution of



1:100 (40x lens).

Immunofluorescence - NFE2L1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using NFE2L1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.