

# PAEP Polyclonal Antibody



Catalog Number: 162248

## Product name

PAEP Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to PAEP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 19-180 of human PAEP (NP\_002562.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - PAEP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using PAEP antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - PAEP Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse pancreas using PAEP antibody at dilution of 1:100 (40x lens).