

PTCH1 Polyclonal Antibody

Catalog Number: 162252

Product name

PTCH1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to PTCH1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 770-970 of human PTCH1 (NP_000255.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IF

Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

Validations




Western blot - PTCH1 Polyclonal Antibody

Western blot analysis of extracts of rat skeletal muscle, using PTCH1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.



Immunofluorescence - PTCH1 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using PTCH1 antibody at dilution of 1:100. Blue: DAPI




for nuclear staining.

PTCH1 Polyclonal Antibody

Catalog Number: 162252

Immunofluorescence - PTCH1 Polyclonal Antibody


Immunofluorescence analysis of L-929 cells using PTCH1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.


Immunofluorescence - PTCH1 Polyclonal Antibody

Immunofluorescence analysis of mouse embryos cells using PTCH1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunofluorescence - PTCH1 Polyclonal Antibody

Immunofluorescence analysis of C6 cells using PTCH1 antibody at dilution of 1:100. Blue: DAPI



for nuclear staining.

Immunofluorescence - PTCH1 Polyclonal Antibody

Immunofluorescence analysis of mouse embryos cells using PTCH1 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.