

SLC6A4 Polyclonal Antibody



Catalog Number: 162256

Product name

SLC6A4 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to SLC6A4

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human SLC6A4 (NP_001036.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - SLC6A4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using SLC6A4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.