

ZNF141 Polyclonal Antibody



Catalog Number: 162262

Product name

ZNF141 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZNF141

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 405-474 of human ZNF141 (NP_003432.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200

IF 1:50 - 1:200

Validations



Immunohistochemistry - ZNF141 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse lung using ZNF141 antibody at dilution of 1:100 (40x lens).