

SLC5A6 Polyclonal Antibody



Catalog Number: 162273

Product name

SLC5A6 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to SLC5A6

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 500 to the C-terminus of human SLC5A6 (NP_066918.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IHC

Dilutions


IHC 1:50 - 1:200

Validations



Immunohistochemistry - SLC5A6 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human liver using SLC5A6 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - SLC5A6 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat heart using SLC5A6 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - SLC5A6 Polyclonal Antibody

SLC5A6 Polyclonal Antibody



Catalog Number: 162273

Immunohistochemistry of paraffin-embedded

mouse heart using SLC5A6 antibody at dilution of 1:100 (40x lens).