

# SLC25A17 Polyclonal Antibody



Catalog Number: 162291

---

## Product name

SLC25A17 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to SLC25A17

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SLC25A17 (NP\_006349.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - SLC25A17 Polyclonal Antibody

Immunofluorescence analysis of 293T cells using SLC25A17 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.