

ZNF562 Polyclonal Antibody



Catalog Number: 162327

Product name

ZNF562 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to ZNF562

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-180 of human ZNF562 (NP_001123503.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


IHC, IF

Dilutions

IHC 1:50 - 1:200


IF 1:50 - 1:200

Validations




Immunohistochemistry - ZNF562 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human lung cancer using ZNF562 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ZNF562 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human gastric cancer using ZNF562 antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - ZNF562 Polyclonal


ZNF562 Polyclonal Antibody



Catalog Number: 162327

Antibody

Immunohistochemistry of paraffin-embedded mouse liver using ZNF562 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZNF562 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using ZNF562 antibody at dilution of 1:100 (40x lens).