

# ZNF415 Polyclonal Antibody



Catalog Number: 162337

## Product name

ZNF415 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to ZNF415

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 100-260 of human ZNF415 (NP\_060825.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - ZNF415 Polyclonal Antibody


Immunohistochemistry of paraffin-embedded rat brain using ZNF415 antibody at dilution of 1:100



(40x lens).

Immunohistochemistry - ZNF415 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse brain using ZNF415 antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - ZNF415 Polyclonal Antibody

# ZNF415 Polyclonal Antibody



Catalog Number: 162337

---

Immunohistochemistry of paraffin-embedded

mouse kidney using ZNF415 antibody at dilution of 1:100 (40x lens).