

# DDX59 Polyclonal Antibody



Catalog Number: 162353

---

## Product name

DDX59 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to DDX59

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human DDX59 (NP\_001026895.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

## Validations



Immunohistochemistry - DDX59 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat testis using DDX59 antibody at dilution of 1:100 (40x lens).