

JAM3 Polyclonal Antibody



Catalog Number: 162357

Product name

Rabbit IgG

JAM3 Polyclonal Antibody

Applications

Specificity

WB

Human

Dilutions


Antibody description

WB 1:500 - 1:2000

Polyclonal antibody to JAM3

Validations

Preparation



Antigen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human JAM3 (NP_116190.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Western blot - JAM3 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using JAM3 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 90s.