

RIMS2 Polyclonal Antibody



Catalog Number: 162387

General Information

Product name:

RIMS2 Polyclonal Antibody

Immunogen:

Applications:

WB

Specificity:

Human, Rat

Antibody description:

Polyclonal antibody to RIMS2

Properties

Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-350 of human RIMS2 (NP_055492.3).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

WB 1:500 - 1:2000

Validations



Western blot - RIMS2 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RIMS2 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 15s.