

# A1BG Polyclonal Antibody



Catalog Number: 162400

## Product name

A1BG Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to A1BG

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 236-495 of human A1BG (NP\_570602.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



## Western blot - A1BG Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using A1BG antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



## Immunofluorescence - A1BG Polyclonal Antibody

Immunofluorescence analysis of 293T cells using A1BG antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.