

ACADSB Polyclonal Antibody



Catalog Number: 162402

Product name

ACADSB Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to ACADSB

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human ACADSB (NP_001600.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - ACADSB Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ACADSB antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.