

# CYP2A6 Polyclonal Antibody



Catalog Number: 162411

---

## Product name

CYP2A6 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to CYP2A6

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 91-300 of human CYP2A6 (NP\_000753.3).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - CYP2A6 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CYP2A6 antibody at dilution of 1:100 (40x lens).  
Blue: DAPI for nuclear staining.