

CYP2J2 Polyclonal Antibody



Catalog Number: 162414

Product name

CYP2J2 Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to CYP2J2

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human CYP2J2 (NP_000766.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

IF

Dilutions

IF 1:50 - 1:200

Validations



Immunofluorescence - CYP2J2 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using CYP2J2 antibody at dilution of 1:100 (40x lens).
Blue: DAPI for nuclear staining.