

SOX10 Polyclonal Antibody



Catalog Number: 162445

General Information

Product name:

SOX10 Polyclonal Antibody

Immunogen:

Applications:

WB

Specificity:

Human

Antibody description:

Polyclonal antibody to SOX10

Properties

Preparation:

Antigen: A synthetic peptide corresponding to a sequence within amino acids 200-300 of human SOX10 (NP_008872.1).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type: Rabbit IgG

Applications

WB 1:500 - 1:2000

Validations



Western blot - SOX10 Polyclonal Antibody

Western blot analysis of extracts of A-375 cells, using SOX10 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.