

# TIMP3 Polyclonal Antibody



Catalog Number: 162446

---

## Product name

TIMP3 Polyclonal Antibody

## Specificity

Human, Mouse, Rat, Monkey

## Antibody description

Polyclonal antibody to TIMP3

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 24-211 of human TIMP3 (NP\_000353.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - TIMP3 Polyclonal Antibody

Immunofluorescence analysis of MCF-7 cells using TIMP3 antibody . Blue: DAPI for nuclear staining.