

EIF5B Polyclonal Antibody



Catalog Number: 162458

Product name

EIF5B Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to EIF5B

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human EIF5B (NP_056988.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - EIF5B Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EIF5B antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.Detection: ECL Basic Kit .Exposure time: 3s.