

COL5A1 Polyclonal Antibody



Catalog Number: 162470

Product name

COL5A1 Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to COL5A1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 20-120 of human COL5A1 (NP_000084.3).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - COL5A1 Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using COL5A1 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.