

EIF4ENIF1 Polyclonal Antibody



Catalog Number: 162491

Product name

EIF4ENIF1 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to EIF4ENIF1

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 150-250 of human EIF4ENIF1 (NP_001157973.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - EIF4ENIF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using EIF4ENIF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 3s.