

MCCC2 Polyclonal Antibody



Catalog Number: 162494

Product name

MCCC2 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to MCCC2

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 284-563 of human MCCC2 (NP_071415.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

Validations



Immunofluorescence - MCCC2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using MCCC2 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.