

# ZNF408 Polyclonal Antibody



Catalog Number: 162496

---

## Product name

ZNF408 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ZNF408

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human ZNF408 (NP\_079017.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000

IF 1:50 - 1:200

## Validations



Immunofluorescence - ZNF408 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ZNF408 antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.