

# UPF1 Polyclonal Antibody



Catalog Number: 162509

## Product name

UPF1 Polyclonal Antibody

## Specificity

Human, Mouse

## Antibody description

Polyclonal antibody to UPF1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 270-370 of human UPF1 (NP\_002902.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:500 - 1:2000


IF 1:50 - 1:200

## Validations



Western blot - UPF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using UPF1 antibody at 1:500 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.



Immunofluorescence - UPF1 Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using UPF1 antibody . Blue: DAPI for nuclear staining.