

# BCLAF1 Polyclonal Antibody



Catalog Number: 162521

## Product name

BCLAF1 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to BCLAF1

## Preparation

Antigen: A synthetic peptide of human BCLAF1

## Formulation

PBS with 0.02% sodium azide, pH7.3.

## Storage

Store at 4°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:1000

IHC 1:50 - 1:100

## Validations



Western blot - BCLAF1 Polyclonal Antibody

Western blot analysis of extracts of MCF7 cells, using BCLAF1 antibody. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST.