

Pan Acetyl-Lysine mAb



Catalog Number: 162531

Product name

Pan Acetyl-Lysine mAb

Specificity

ALL

Antibody description

Monoclonal antibody to Pan Acetyl-Lysine

Preparation

Antigen: Acetylated BSA

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Monoclonal

Ig Type

Mouse IgG

Applications

WB

Dilutions

WB 1:500 - 1:1000

Validations



Western blot - Pan Acetyl-Lysine mAb

Western blot analysis of extracts of HeLa cells, using Pan Acetyllysine antibody .