

ATP1A4 Polyclonal Antibody



Catalog Number: 162535

Product name

Rabbit IgG

ATP1A4 Polyclonal Antibody

Applications

Specificity

WB

Mouse, Rat

Dilutions

Antibody description

WB 1:200 - 1:2000

Polyclonal antibody to ATP1A4

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-90 of human ATP1A4 (NP_653300.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Validations



Western blot - ATP1A4 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using ATP1A4 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1S.