

CEBPD Polyclonal Antibody



Catalog Number: 162538

Product name

CEBPD Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to CEBPD

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-269 of human CEBPD (NP_005186.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - CEBPD Polyclonal Antibody

Western blot analysis of extracts of mouse lung, using CEBPD antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90S.