

TMEM55A Polyclonal Antibody



Catalog Number: 162663

Product name

TMEM55A Polyclonal Antibody

Specificity

Mouse

Antibody description

Polyclonal antibody to TMEM55A

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 130-195 of human TMEM55A (NP_061180.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - TMEM55A Polyclonal Antibody

Western blot analysis of extracts of mouse liver, using TMEM55A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 5s.