

# VMP1 Polyclonal Antibody



Catalog Number: 162684

## Product name

VMP1 Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to VMP1

## Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 300 to the C-terminus of human VMP1 (NP\_112200.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IF

## Dilutions

WB 1:200 - 1:2000

IF 1:50 - 1:200

## Validations



Western blot - VMP1 Polyclonal Antibody

Western blot analysis of extracts of mouse kidney, using VMP1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 5s.