

OTUD7A Polyclonal Antibody



Catalog Number: 162717

Product name

OTUD7A Polyclonal Antibody

Specificity

Human

Antibody description

Polyclonal antibody to OTUD7A

Preparation

Antigen: A synthetic peptide corresponding to a sequence within amino acids 800-900 of mouse OTUD7A (NP_570950.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:200 - 1:2000

Validations



Western blot - OTUD7A Polyclonal Antibody

Western blot analysis of extracts of U-251MG cells, using OTUD7A antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.