

# Acetyl-Histone H2A-K5 pAb



Catalog Number: 162737

## Product name

Acetyl-Histone H2A-K5 pAb

WB 1:500 - 1:2000

IHC 1:50 - 1:100

IF 1:50 - 1:100

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to Acetyl-Histone H2A-K5

## Preparation

Antigen: A synthetic acetylated peptide around K5 of human Histone H2A (NP\_003508.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications

WB, IHC, IF


## Dilutions

## Validations



Immunohistochemistry - Acetyl-Histone H2A-K5 pAb

Immunohistochemistry of paraffin-embedded human breast using H2AK5ac antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - Acetyl-Histone H2A-K5 pAb

Immunohistochemistry of paraffin-embedded mouse liver using H2AK5ac antibody at dilution of 1:100 (40x lens).