

BMP6 Polyclonal Antibody



Catalog Number: 162748

General Information

Product name:

BMP6 Polyclonal Antibody

Immunogen:

Applications:

IF

Specificity:

Human

Antibody description:

Polyclonal antibody to BMP6

Properties

Preparation:

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 375-513 of human BMP6 (NP_001709.1).

Formulation:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage:

Store at -20°C. Avoid freeze / thaw cycles.

Clonality:


Polyclonal

Ig Type:Rabbit IgG

Applications

IF 1:50 - 1:200

Validations



Immunofluorescence - BMP6 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using BMP6 antibody at dilution of 1:100. Blue: DAPI for nuclear staining.